0476

RI.MOS. s.r.l. viale Gramsci, 29 – 41037 MIRANDOLA (MO) ITALY – tel. +39 0535 25755 –
fax +39 0535 25509 - e-mail: info@rimos.com – www.rimos.com

GINRAM®

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E

Indications: this device is intended for the control or reduction of postpartum uterine hemorrhage in cases where conservative therapeutic management is indicated.

INSTRUCTIONS FOR USE:

Before the transvaginal or transabdominal positioning of the postpartum Metrastop balloon, remove all placental fragments from the uterus and examine the patient in order to exclude any laceration or trauma in the genital tract, then ensure that the bleeding is not of arterial origin.

Transvaginal placement after spontaneous delivery

1-determine the uterine volume by direct examination or ultrasound examination.

2-Insert the distal part of the catheter into the uterus, making sure that the balloon is completely positioned beyond the cervical canal and the internal ostium.

It is recommended to use ring forceps to position on the uterine cervix in order to perform a slight traction and favor the introduction of the device under vision.

If this is not possible, it is advisable to grab the device near the lumen-balloon seam and accompany its introduction into the cervix by hand.

3-insert an urinary catheter in order to collect and monitor urine production.

The control of correct positioning of the device must be under ultrasound vision.

Transabdominal placement during C-section

1-determine the uterine volume by direct examination.

2-Insert the catheter by hysterotomy, placing the part with inflation fittings into the vaginal canal through cervix, and the tip with balloon into the uterine cavity.

Carry out, by vaginal pathway, a control of the access to the inflation lumens.

3- Check the correct balloon positioning inside the uterine cavity.

4- Close the incision as standard practice, ensuring to avoid the sting of the balloon during the surgical suture.

5-Check of the suture keeping and evolution of the bleeding during the device inflation under direct vision of the professional.

NOTE- before the balloons inflation, ensure that all the components of the product are intact and hysterotomy has been properly sutured.

If clinically relevant, abdomen can remain open while inflating the balloons, so as to monitor carefully uterine distension and to confirm the closure of the hysterotomy.

NOTE- if clinically relevant, a B-Lynch compression suture can be carried out in conjunction with Metrastop postpartum device.

Transvaginal placement during C-section

In the event of any hemorrhage after the hysterotomy suture, it is possible to position Metrastop device through the vaginal canal and uterine cervix monitoring the manoeuvre under direct vision of the professional, avoiding the laceration of hysterotomic suture.

Verify the suture keeping by direct examination after the balloons inflation.

Balloons inflation – through 60 ml syringe

The external balloon has been tested up to 60 ml saline solution inflation.

The internal balloon has been tested up to max. 420 cc air inflation (7 pumps).

Reccomended inflation: 350 cc of air (6 pumps).

1-Preload the 60 ml included syringe with sterile saline solution and connect it to the blue lumen named “water”.

Inject the content of the syringe and close the stopcock.

In this way the external balloon is inflated.

2-Connect the same 60 ml syringe (empty) to the inflation valve of the stopcock of the transparent lumen named “air inflation”, then insufflate max 420 cc air (7 pumps).

Reccomended inflation: 350 cc of air (6 pumps).

It is advised to repeat loud the number of times it’s done the air inflation manoeuvre (e.g. 1, 2 etc. until 7).

Close the stopcock after the inflation.

The valve linked to the transparent lumen “air inflation” enables to keep the syringe connected during the inflation procedure without unplugging it.

In this way the internal balloon is inflated.

3-Connect the red lumen for the drainage to the included fluid collection bag in order to monitor the hemostasis.

In order to do an appropriate monitoring of hemostasis, the balloon drainage connector and the lumen can be released from the blood clots by wash with sterile saline solution.

In case of complete filling of the collection bag, it is possible to replace it with an other one supplied by th ward (e.g. collection bag for bladder catheter).

After 4 hours since the positioning, it is advised an ultrasound evaluation of the intra-uterine balloon placement (both internal and external), as well as the verification of the established volume maintenance during placement.

In case of volume loss of one of the two balloons, it is advised to use the blue lumen (water) in order to restore the volume needed by injection of sterile saline solution.

Device removal

The device can be removed in advance after that the physician determines the obtaining of haemostasis .

Metrastop must be kept in place no longer than 24 hours.

1-Remove the inflation valve of the transparent lumen and open the stopcock in order to leak the air out. In this way the internal balloon is deflated.

2-Open the blue lumen stopcock or use a suitable syringe to draw the contents of the external balloon until it totally deflates.

It is also possible to cut the catheter lumens in order to deflate the balloons more quickly.

**0476**

RI.MOS. s.r.l. viale Gramsci, 29 – 41037 MIRANDOLA (MO) ITALY – tel. +39 0535 25755 – fax +39 0535 25509 – e-mail: info@rimos.com – www.rimos.com

3-It's advised to wait some minutes in order to verify the effective absence of any pathological uterine bleeding before to utterly remove the device.

Subsequently pull softly the catheter out of the uterus and vaginal canal, then throw it.

4-monitor the patient in order to detect potential signs of bleeding resuming.

Contraindications

Arterial bleeding requiring surgical exploration or angiographic embolization.

Cases indicating a hysterectomy

Pregnancies

Cervical cancer

Vaginal, cervix or uterine purulent infection

Untreated uterine anomaly

Disseminated intravascular coagulation

Warnings

Read the instructions before use.

Continue to monitor the patient after Metrastop application through clinical, laboratory and instrumental assessment in order to identify signs of any persistent uterine bleeding and/or other sources of anemization.

Even after Metrastop removal, control any potential bleeding persistence.

Metrastop is not a substitute for surgical treatment, whenever indicated.

Metrastop is not a substitute for any therapy appropriate to achieve the blood volume (fluids, plasma, blood).

After 4 hours since the positioning, it is advised an ultrasound evaluation of the intra-uterine balloon placement (both internal and external), as well as the verification of the established volume maintenance during placement.

In case of volume loss of one of the two balloons, it is advised to use the blue lumen (water) in order to restore the volume needed by injection of sterile saline solution.

There are no clinical data on Metrastop use in disseminated intravascular coagulation (DIC). During Metrastop usage it is necessary to monitor the urine output through the use of a bladder catheter.

Assess the clinical utility of Metrastop device under conditions of postpartum uterine hemorrhage determined by altering placentation (placenta accreta, increta, percreta).

Do not reuse. Use only if the packaging is intact. Use immediately after the pouch opening.

This device must be used only by specialist medical staff.

Do not use the device in the following cases: DIC, ongoing infections, arterial bleeding that needs surgical therapy.

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LOT BATCH CODE	REF REFERENCE NUMBER	ATTENTION. SEE INSTRUCTIONS FOR USE	LATEX FREE	STERILE EO STERILIZED BY ETHYLENE OXIDE
USE BY MM/YYYY	CONSULT INSTRUCTIONS FOR USE	MANUFACTURER	DATE OF MANUFACTURE	DO NOT REUSE