



0476



RI.MOS. s.r.l. viale Gramsci, 29 – 41037 MIRANDOLA (MO) ITALY – tel. +39 0535 25755 –
fax +39 0535 25509 - e-mail: info@rimos.com – www.rimos.com

GINRAM®

PRODOTTI MONOUSO STERILI PER LA GINECOLOGIA

Gynaecological disposable sterile devices

U-CUT

Forbice lunga in plastica con lama in acciaio

Long plastic scissors with steel blades

COD. 770114


I

Istruzioni per l'uso

Impugnare la forbice e utilizzarla per la rimozione di suture, procedure mediche generiche e ove è necessario recidere fili, garze e affini.

Le lame in acciaio inox accuratamente affilate permettono un taglio preciso e sicuro.

Avvertenza: Leggere le istruzioni prima dell'uso. Non riutilizzare. Utilizzare solo se la confezione è integra. Utilizzare immediatamente dopo l'apertura della busta. Distruggere dopo l'utilizzo in accordo alle norme vigenti. Il dispositivo deve essere utilizzato solo da personale medico specializzato.



GB

Instructions for use

Hold the scissor and use it for suture removal, general medical procedures and where it is necessary to cut strings, gauzes and similar.

Stainless steel blades with accurate sharpening allow a precise and safe cut.

Warning: Read the instructions before using the device. Do not reuse. Use only if the package is not damaged. Use immediately after opening. Dispose after use according to the laws in force. For use only under trained medical staff.